



Airflow²

Educational Ducted

This Educational Ducted Fume Cupboard has the capability of removing 100% of the fume and smoke created by chemical reactions.

This Fume Cupboard offers the best protection when undertaking wet chemistry experiments which is ergonomically designed for the capture and containment of fume and smoke.

Available in 1000mm 1200mm and 1500mm widths with a chemically resistant Trespa base and sliding glass sash this is one of the best on the market.

Standard height 1880mm

Compliant with: - Building Bulletin 88 (Rev G9) Designed and manufactured in accordance with Energy efficient design and operation of Fume Cupboards Guide 320, Meets BREAMM Hea 17 and Ene 11 compliance.

Please note a suitable Extraction System will be required for the Fume Cupboard to operate correctly, the Extract System should be designed and installed to meet the requirements of D.W. 154 and Building Bulletin 88.

This Fume Cupboard is 100% designed and built in the UK.



The Fume Cupboard Company

15 Leywood Close, Braintree, Essex CM7 3NP

info@thefumecupboardcompany.co.uk | www.thefumecupboardcompany.co.uk

Telephone: 01376 619371



Airflow 2

Educational Ducted

Options Available:

- Docking Station
- Cast epoxy resin work surface
- Fixed position (no castors)

Standard Features:

- Robust powder coated mild steel construction.
- Heavy duty castors (two front castors braked)
- Trespa Toplab BASE bonded work base 6mm toughened glass rear and side panels.
- 6mm clear acrylic front hinged visor
- Circular slimline LED low energy light fitting.
- Water and Gas hoses are all leak tested by the hose manufacturer prior to delivery to site
- All brass quick release connections are genuine Snap-Tite fittings
- Standard Services & Controls
- Cold water tap and front control valve with flexible hose (2 metres)
- 100mm dia Vulcathene drip cup c/w bottle trap with flexible drain
- Gas tap and front control valve with flexible hose (2 metres)
- Double switched 13A socket outlet
- Light switch
- Fan switch
- Hour Counter
- RCD protection
- Low Airflow Alarm Panel

Specification:

- G9 Revision of DfEE Building Bulletin 88
- Fume cupboard commissioned pre-tested in the factory prior to delivery.



The Fume Cupboard Company

15 Leywood Close, Braintree, Essex CM7 3NP

info@thefumecupboardcompany.co.uk | www.thefumecupboardcompany.co.uk

Telephone: 01376 619371