

Introducing the Advanced Filtered Series

Industrial fume hoods in 7 sizes





A cutting-edge line of ductless fume cupboards crafted for optimal safety and environmental protection.

Choose from our range of 7 standard sizes.

Available in both metal and polypropylene construction.

These high-efficiency solutions are specifically engineered to safeguard users and the surrounding environment from potentially harmful vapours generated during work processes.

Elevate your workspace with the Advanced Filtered Series for unparalleled performance and peace of mind.



Key Features

- High operator protection to fume and particle hazards
- Improved clamping eliminates bypass leakage
- Filter blockage alarm
- Polypropylene work surface (stainless steel option)
- High-capacity filters for more demanding applications

Dimensions

- All models – 1346mm height, 695mm deep
- Heights – 762mm, 864mm, 991mm, 1245mm, 1499mm, 1753mm, 2438mm